

SAS™ Legionella Test

The SAS™ Legionella Test is a visually read, in vitro immunochromatographic rapid assay for the presumptive qualitative detection of *Legionella pneumophila* serogroup 1 antigens in human urine.

For *In-Vitro* Diagnostic Use

Store at 15° to 30° C

For Technical Assistance Call 800-272-2710
Outside the USA Call 210-699-8800



SA Scientific, Ltd.

4919 Golden Quail, San Antonio, Texas 78240

www.sascientific.com

INTENDED USE

The SAS™ Legionella Test is a visually read, in vitro immunochromatographic rapid assay for the presumptive qualitative detection of *Legionella pneumophila* serogroup 1 antigens in human urine. This test is intended to aid in the presumptive diagnosis of Legionnaires' disease in conjunction with culture and other methods for patients with signs and symptoms of pneumonia. This test is for prescription use only.

INTRODUCTION

Legionella pneumophila serogroup 1 is gram-negative bacillus that appears off-white in color and circular in shape. This bacterium is associated with Pontiac fever and Legionnaires' disease. An estimated 25,000² to 10,000³ cases involving *Legionella pneumophila* infections occur in the United States annually. About 5% to 30% of people with Legionnaires' disease die.

Infections caused by this bacterium include the pulmonary disease of pneumonia. In addition, extra-pulmonary infections may be spread through the bloodstream or the lymph system, infecting the heart, brain, kidney, liver, and/or spleen. Early onset symptoms include fever, chills, muscle aches, and lack of appetite. Respiratory symptoms such as coughing or congestion are usually absent. As the disease progresses, a dry, hacking cough develops and may become productive after a few days. Half of people who develop Legionnaires' disease suffer shortness of breath. In many cases¹, a fever as high as 104° F has been reported.

Common methods for the detection of *Legionella pneumophila* are culture, direct fluorescent antibody (DFA), indirect fluorescent antibody (IFA)⁵, enzyme immunoassays (EIA)⁵ and polymerase chain reaction (PCR)^{6,8} using respiratory specimens of expectorated sputum, bronchial washing, transtracheal aspirate and lung biopsy. *Legionella pneumophila* serogroup 1 antigen has been detected in urine during the acute phase of the disease and has become the method by choice for laboratory diagnosis.^{4,5} The SAS™ Legionella test is designed for the detection of *Legionella pneumophila* serogroup 1 soluble antigen in urine.

PRINCIPLE OF THE TEST

The SAS™ Legionella test utilizes a combination of polyclonal antibodies against the antigens of *Legionella pneumophila*. The SAS™ Legionella test begins with the addition of urine to the test device. The specimen is absorbed by the sample pad and then moves through the conjugate pad which contains dried gold conjugated antibodies which are specific for *Legionella pneumophila* antigens; if the Legionella antigens are present in the urine sample, a "half-sandwich" immunocomplex is formed. This immuno-complex then migrates via capillary action along a nitrocellulose membrane containing immobilized antibodies to *Legionella pneumophila* antigens. The immobilized antibodies bind the "half-sandwich" immuno-complex to form a "whole sandwich" immuno-complex. Thus, when the "whole sandwich" is formed, a visible, pink colored line develops in the specimen zone on the test device. In the absence of a Legionella antigen, a "sandwich" immuno-complex is not formed and a negative result

is indicated. To serve as a procedural control, a pink colored control line will always appear in the control zone regardless of the presence or absence of Legionella antigen. The test is available in cassette format.

MATERIALS PROVIDED

1. Test Devices – A membrane coated with a combination of polyclonal antibodies against the antigens of *Legionella pneumophila* serogroup 1.
2. Disposable specimen dropper

MATERIALS NOT PROVIDED

1. Timer
2. Specimen collection containers

PRECAUTIONS

1. For *In-vitro* diagnostic use only.
2. Use pipette provided with kit only.
3. Use separate pipettes for each sample or control.
4. The test devices should remain in the pouch until ready for use.
5. All specimens should be considered potentially hazardous, and should be handled in the same manner as an infectious agent.
6. Wear disposable gloves while handling samples and wash hands after the assay is complete. Warning: The user should refer to the relevant section of the CDC-NIH manual "Biosafety in Microbiology and Biomedical Laboratories," Fifth Edition, 2007.
7. Avoid any contact with eyes.
8. The test device and all materials should be discarded in a proper biohazard container after testing.
9. Do not use kits after the stated expiration date.
10. Do not reuse test devices or kit materials.

STORAGE CONDITIONS

SAS™ Legionella test devices should be kept at room temperature (15°-30°C) in sealed pouches with a desiccant. Do not freeze the test kit.

SPECIMEN COLLECTION AND PREPARATION


The urine specimen must be collected in a clean, sterile container. Urine specimens may be refrigerated (2°-8°C) and stored up to 72 hours prior to assay. Specimens may be stored frozen for up to 2 years. If specimens are refrigerated, they must be equilibrated to room temperature before testing.

TEST PROCEDURE

Allow specimen and/or controls to reach room temperature (15° – 30°C) prior to testing.

1. Remove test device from pouch and lay on a flat surface. Label the device with the patient ID.

2. Dispense 4 drops (approximately 0.15 ml) of specimen into the round sample well (see illustration below). Wait for colored lines to appear.




3. Read results at 15 minutes or less. The presence of the control line is not indicative of the test being completed. Interpretation after 15 minutes may be inaccurate.

INTERPRETATION OF RESULTS


Negative Results

The test is negative if a colored line appears in the C area (control) only.



Positive Results

The test is positive if one colored line appears in the S (specimen) area and one colored line appears in the C (control) area. Any colored line in the specimen area should be considered positive. Colored lines may appear lighter or darker than one another. Specimens with Legionella antigen levels near the threshold of the test may develop color (faint lines) overtime after the 15 minute reading; however, a line that appears after 15 minutes is not diagnostic and should be ignored.



Invalid Results

The test is invalid if no colored line appears in the C (control) area even if a colored line appears in the S (specimen) area. If a colored line does not appear in the C area, the test is invalid and

should be repeated. Colored lines that appears in the S (specimen) area after 15 minutes are not diagnostic and should be ignored.

Quality Control Internal Controls

The appearance of a Control line in the C region of the device is a positive procedural control. Correct procedural technique, specimen flow and device performance is confirmed when a colored line appears in the C (control) area of the membrane. If the colored line fails to appear in the C (control) area, the test result is invalid.

A clear background is an internal negative procedural control. The background color should be white to light pink and should not interfere with the reading of the test result. If any result is difficult to interpret, the test should be repeated with the same sample to eliminate the potential for error. Obtain a new sample and retest when the original sample repeatedly produces unreadable results. It is recommended that the internal quality controls be recorded for each sample run.

External Controls

Quality control should be performed on each new test kit box opened. Additional testing of the appropriate controls of *Legionella pneumophila* may be performed in accordance with federal, state and/or local regulations, accrediting requirements, and your laboratory's quality control procedures. Positive and negative external controls are supplied separately.

LIMITATIONS OF THE PROCEDURE

- A negative antigen result does not rule out the possible infection from *Legionella pneumophila* serogroup 1. The antigen concentration may fall below the detectable limit of the test. Culture is recommended for patients with suspected pneumonia to determine if other causative agents other than *Legionella pneumophila* serogroup 1 are responsible.
- The test provides a presumptive diagnosis for only *Legionella pneumophila* serogroup 1 infection. In order to make an accurate diagnosis, culture results, serology testing and antigen detection methods in conjunction with clinical findings should be used.
- Sensitive immunoassays may demonstrate positive results with specimens containing heterophilic antibodies. If the qualitative interpretation is inconsistent with clinical findings, then further testing by an alternate method should be performed.
- Excretion of Legionella antigen may be seen as early as three (3) days after the onset of symptoms. However, antigen may continue to be excreted up to one (1) year afterwards⁷. A positive result with the SASTM Legionella Test may occur during a current and past infection; therefore, in order to make an accurate diagnosis, culture results, serology testing and antigen detection methods in conjunction with clinical findings should be used.

PERFORMANCE DATA Clinical Sensitivity and Specificity (Retrospective Study)

Three clinical sites (USA and Netherlands) tested three hundred twenty four (324) frozen specimens using the SASTM Legionella Cassette test. These samples were previously tested for Legionella by cell culture. Results are reported below.

SAS TM Legionella Test		Cell Culture		
		+	-	
	+	95	11	106
	-	10	208	218
		105	219	324

Sensitivity: 95/105 x 100 = 90.5 % (95% CI: 83.4-94.8)
Specificity: 208/219 x 100 = 95.0 % (95% CI: 91.2-97.5)

Analytical Specificity

Forty-Nine (49) fresh patient urine samples from healthy individuals were collected and assayed at one (1) clinical site. One hundred percent (100%) of these were found to be negative by the SASTM Legionella test.

Cross-Reactivity with other Respiratory Tract Infections

Ninety-nine (99) urines from patients diagnosed for other etiological respiratory tract infections (84 culture confirmed, 15 suspected) were tested using the SASTM Legionella Test. The results showed a lack of reactivity in 98/99 samples (99.0%).

Bacterial Cross-Reactivity and Interference

To confirm the analytical specificity of the SASTM Legionella Test, bacterial cultures likely to be found in the respiratory tract were tested. Bacterial cultures were tested at 1 x 10⁸ cfu/ml. All yielded negative results.

To confirm a lack of interference by other bacterial species in the SASTM Legionella Test, purified Legionella antigen was added to bacterial cultures likely to be found in the respiratory tract. The concentration of Legionella antigen was 1 x 10⁵ cfu/ml. All tests yielded positive results.

<i>Legionella pneumophila</i> SG2	<i>Legionella gormanii</i>	<i>Proteus mirabilis</i>
<i>Legionella pneumophila</i> SG3	<i>Legionella longbeachae</i>	<i>Citrobacter freundii</i>
<i>Legionella pneumophila</i> SG4	<i>Legionella feeleii</i>	<i>Enterobacter cloacae</i>
<i>Legionella pneumophila</i> SG5	<i>Staphylococcus aureus</i>	<i>Serratia marcescens</i>
<i>Legionella pneumophila</i> SG6	<i>Staphylococcus epidermidis</i>	<i>Pseudomonas aeruginosa</i>
<i>Legionella bozemanii</i>	<i>Enterococcus faecalis</i>	<i>Candida albicans</i>
<i>Legionella micdadei</i>	<i>Escherichia coli</i>	<i>Streptococcus pneumoniae</i>
<i>Legionella dumoffii</i>	<i>Klebsiella pneumoniae</i>	<i>Haemophilus influenzae</i>

Common Urine Components and Drugs

The components and drugs listed below were tested to determine if they interfered or cross-reacted with the SASTM Legionella Test. Extracted *Legionella pneumophila* antigen was added at a concentration of 4.5 x 10⁴ cfu/mL. In addition, each compound was tested in the absence of *Legionella pneumophila* antigen. None listed interfered or cross-reacted with the test results.

Glucose	2000mg/dL	Protein - BSA	500mg/dL
Urea	2000mg/dL	Bilirubin	20mg/dL
Ascorbic Acid	100mg/dL	Rifampicin	0.09mg/mL
Erythrocytes	10 ⁶ /mL	Erythromycin	0.067mg/mL
Leukocytes	10 ⁶ /mL	Ciprofloxacin	0.22mg/mL

Limit of Detection

The limit of detection of the SASTM Legionella test was determined to be 5 x 10⁴. *Legionella pneumophila* serogroup 1 ATCC 33152 was prepared using BCYE agar. A dilution of the working concentration was performed. The limit of detection of the SASTM Legionella Cassette Test was determined from these concentrations. The test was performed in 2 replicates. Results were recorded at 15 minutes.

Reproducibility Study










The reproducibility of the SASTM Legionella Test was evaluated at three clinical laboratory sites. The SASTM Legionella Test was tested against a panel of six (6) specimens of which included four levels of positives and two negatives. The low and high positives were from the purified Legionella antigen. Negatives were comprised of either urine or Legionella antigen below the detectable limit. Three (3) different laboratory personnel assayed each specimen at each laboratory facility over 3 day s. The overall reproducibility for the SASTM Legionella Test was 100%.

REFERENCES

- Edelstein, P.H. and R.D. Meyer. "Legionella Pneumonias." In Respiratory Infections: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
- Marston, B.J., H. B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch. Intern. Med 1994;154:2417:2422
- Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Breiman. Prospects for vaccine development. Presented at the 4th International Symposium on *Legionella*, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. *Legionella: Current status and emerging perspectives*. Washington, DC American Society for Microbiology, 1993.
- Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
- Lindsay, D. S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S Edwards. Laboratory diagnosis of legionnaires' disease due to *Legionella pneumophila* serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.

- Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of *Legionella* species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Microbiol. 38; 2000; 1709-1712.
- Kohler, R.B., W.C. Winn, L.J. Wheat. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
- Jonas, Daniel, A. Rosenbaum, S. Weirick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Microbiol. 33; 1995; 1247-1252.

KEY TO INTERNATIONAL SYMBOLS in ENGLISH

SYMBOLS	FOR
	Manufacturer
	Expiration Date
	In-vitro diagnostic device
	Authorized Representative in European Community
	Contains sufficient for <n> tests
	Batch Code
	Catalog Number
	Product conforms to IVD directive
	Temperature Limitation (15-30°C)

Authorised Representative:
MegaCor GmbH
Europaplatz 1
88131 Lindau
Germany



Rev 09-14
70-PI-LEGIONELLA