

# Quick Reference Guide

Description	Cat. No.	Test/ Kit	Specimen	Application/Shelf Life	Sensitivity, Specificity	CLIA/ CPT Code	Results/ Minutes
-------------	----------	--------------	----------	------------------------	-----------------------------	-------------------	---------------------

## FERTILITY

Human Chorionic Gonadotropin (hCG)							
SAS™ Pregnancy Strip	087525	25	Urine	Detects hCG in urine (Shelf Life: 24 mos.)	25 mIU/mL, >99%	Waived/81025	4
Urine Control Set	087413	25 mIU, 250 mIU, Negative, 3 Vials/Pk (NOTE: For use with 087525, 087430, 087730, 087750, 087731, and 087751)					
SAS™ One-Step Pregnancy	087430	30	Urine	Detects hCG in urine (Shelf Life: 24 mos.)	25 mIU/mL, >99%	Waived/81025	4
Urine Control Set	087413	25 mIU, 250 mIU, Negative, 3 Vials/Pk (NOTE: For use with 087525, 087430, 087730, 087750, 087731, and 087751)					
SAS™ Serum/Urine hCG	087730	30	Serum/Urine	Detects hCG in serum or urine (Shelf Life: 21 mos.)	25 mIU/mL, >99%	Serum: Non-waived/84703 Urine: Waived/81025	4
	087750	50					
Serum Control Set	087712	25 mIU, 250 mIU, Negative, 3 Vials/Pk (NOTE: For use with 087730, 087750, 087731, and 087751)					
SAS™ Ultra hCG Serum/Urine	087731	30	Serum/Urine	Detects hCG in serum or urine (Shelf Life: 18 mos.)	Serum: 10 mIU/mL & Urine: 20 mIU/mL, >99%	Serum: Non-waived/84703 Urine: Waived/81025	Serum: 5-8, Urine: 4
	087751	50					
Serum Control Set	087715	10 mIU, 250 mIU, Negative, 3 Vials/Pk (NOTE: For use with 087731, and 087751)					

## INFECTIOUS DISEASES

Viral							
SAS™ Adeno Test	067020	20	Cell Culture/Eye Swab/ Nasopharyngeal Secretions/ Fecal Material	Detects all Adenovirus strains (Shelf Life: 18 mos.)	96.8%, 94.1%	Non-waived/87809	15
Control Set	067101	2 x 1 mL, 2 Vials/Pk (NOTE: For use with 067020)					
SAS™ Rota Test	069020	20	Stool	Detects wide range of Rotavirus serotypes (Shelf Life: 18 mos.)	>98%, >98%	Non-waived/87899	15
Control Set	069101	2 x 1 mL, 2 Vials/Pk (NOTE: For use with 069020)					
SAS™ RSVAlert (Respiratory Syncytial Virus)	041111	10	Nasopharyngeal Washes/ Aspirates & Swabs	Detects RSV antigen, strain types A & B (Shelf Life: 18 mos.)	97%, 98%	Waived/87807QW* Non-waived/87807	15
	041130	30					
Control Set	041101	2 x 1 mL, 2 Vials/Pk (NOTE: For use with 041111 and 041130)					
SAS™ FLUAlert (30 Test Combo)	046230	15/Ea. A & B	Nasopharyngeal Washes & Aspirates	Detects Influenza A & B (Shelf Life: 18 mos.)	Flu A: 76%, 97.9%, Flu B: 90.5%, >99%	Waived/87804*†	15/Each A & B
Control Set	046213	3 x 1 mL, 1 Vial Each of Flu A, Flu B, and Negative (NOTE: For use with 046230)					

BACTERIAL							
SAS™ StrepAlert® (Dipstick)	611530	30	Throat Swab	Detection of Strep A (Shelf Life: 18 mos.)	97.9%, 97.3%	Waived & Non-waived/ 87880QW* & 87880	5
SAS™ Legionella	044020	20	Urine	Detects <i>Legionella pneumophila</i> serogroup 1 antigens (18 mo. shelf life)	> 90.5% >95%	Moderately complex/ 87499	15
Control Set	044413	2 x 1 mL, 2 Vials/Pk (NOTE: For use with 044020)					

## SICKLE CELL

SAS™ Sickle Cell	094100	100	Whole Blood	Detects Homozygous (S/S) & Heterozygous (A/S) Sickle Cell (Shelf Life: 24 mos.)	>99%, >99%	Non-waived/85660	5
------------------	--------	-----	-------------	---	------------	------------------	---

\* "QW" modifier is added to report use of CLIA-waived test systems for Medicare/Medicaid claims.

† The SAS™ Flu A and B are run separately. Consequently, reimbursement can be claimed when running each Flu A and Flu B assay.

Distributed by:

Manufactured by:



Reliable Diagnostic Tests

4919 Golden Quail, San Antonio, Texas 78240 U.S.A.

(210) 699-8800 • (210) 699-6545 Fax • www.sascientific.com